

NEPHROSTOMY CATHETER SET (INFANT PCN SET)



Set Includes

Initial Puncture Needle,
Guidewire, Dilators,
Pigtail Catheter,
Scalpel Blade &
Urine Bag Connector

Size (in Fr)

6

Intended Use:

Used for percutaneous placement of a pigtail catheter in the renal Pelvis for nephrostomy drainage in neonates and infants.

Key Features:

- Designed for PCN in infants with small pelvis where the standard Pigtail catheters are too big to place
- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- Radiopaque for radiologic confirmation of the position
- Biocompatible tested material
- Measure marking on the needle and the catheter