

OPTICORE[®]

BIOPSY GUN - FULLY AUTOMATIC

Used for obtaining histological core samples from soft tissue and organs.

Key Features

- Echogenic and atraumatic needle tip
- Ergonomic and light weight
- Cannula with depth marking
- Easy and one hand operated with two stage arming
- Available with/without Coaxial Needle



Biopsy Gun - Fully Automatic

<i>Size</i>	<i>Length</i>
14 G	10 & 16 cm
16 G	10, 16 & 20 cm
18 G	10, 16, 20 & 25 cm
20 G	10, 16 & 20 cm