


PCN Catheter With Locking Mechanism (PCN Catheter With Suture)



Package Includes
PCN Catheter with locking mechanism, metal or plastic straightener & urine bag connector

Size (in Fr)
6, 7, 8, 10, 12, 14 & 16

Length (in cm)
20, 22, 25, 30, 35 & 40

Intended Use:

PCN Catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine. The locking mechanism prevents accidental slippage of the catheter.

Key Features:

- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- The locking mechanism prevents accidental slippage of the catheter
- Eye holes for better drainage
- Also available in Hydrophilic coating for easy insertion

Note: Once the catheter is deployed pull the string, wind and lock to the hub of the catheter.