

## PCN Catheter With Needle



**Package Includes**  
PCN Catheter with metal straightener, metal trocar & urine bag connector

### Intended Use:

Used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage by direct puncture.

### Key Features:

- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- Eye holes for better drainage
- Trocar design and catheter design prevents buckling of catheter during the procedure
- Also available in Hydrophilic coating for easy insertion
- Also available with locking mechanism

### PCN Catheter with Needle

#### Size (in Fr)

6, 7, 8, 8.5, 9, 10, 12, 14, 16 & 18

#### Length (in cm)

20, 22, 25, 30, 35 & 40

### PCN Catheter with Needle & Locking Mechanism

#### Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

#### Length (in cm)

20, 22, 25, 30, 35 & 40