


PCN Set - Malecot (Malecot Catheter Set)



Set Contains

- IP Needle
- Guidewire J tip
PTFE coated
- Fascial Dilators
- Urine bag connector
- Malecot Catheter
- Scalpel blade
- Catheter fixing disc

Size (in Fr)

8, 10, 12, 14, 16, 18, 20 & 22

Length (in cm)

20, 25, 30 & 35

Intended Use:

Used for percutaneous placement of a Malecot Catheter in the renal pelvis for nephrostomy drainage.

Advantage:

Kit provides entire devices for performing Nephrostomy procedure.