


## PCN Set - Pigtail With Locking Mechanism (Nephrostomy Catheter Set)



### Set Contains

- IP Needle
- Guidewire J tip  
PTFE coated
- Fascial Dilators
- PCN Catheter with  
Locking Mechanism
- Urine bag connector
- Scalpel blade
- Catheter fixing disc

### Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

### Length (in cm)

20, 22, 25, 30, 35 & 40

### Intended Use:

Used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage.

### Advantage:

Kit provides entire devices for performing Nephrostomy procedure