

# PCN CATHETER

PCN catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine.

## Key Features

- The pigtail tip configuration helps in the retention of the distal tip in renal pelvis
- Eye holes for better drainage
- Also available with Hydrophilic coating for easy insertion

## Package Includes

PCN Catheter, metal or plastic straightener & urine bag connector

## Size (in Fr)

6, 7, 8, 8.5, 9, 10, 12, 14, 16 & 18

## Length (in cm)

20, 22, 25, 30, 35 & 40

