

PCN CATHETER WITH LOCKING MECHANISM

PCN catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine. The locking mechanism prevents accidental slippage of the catheter.

Key Features

- The pigtail tip configuration helps in the retention of the distal tip in renal pelvis
- The locking mechanism prevents accidental slippage of the catheter
- Eye holes for better drainage
- Also available with Hydrophilic coating for easy insertion

Package Includes

PCN Catheter with locking mechanism, metal or plastic straightener & urine bag connector

Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

Length (in cm)

20, 22, 25, 30, 35 & 40



Note: Once the catheter is deployed, pull the string, wind and lock to the hub of the catheter.