

# PCN CATHETER WITH NEEDLE

Used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage by direct puncture.

## Key Features

- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- Eye holes for better drainage
- Trocar design and catheter design prevents buckling of catheter during the procedure
- Also available in Hydrophilic coating for easy insertion
- Also available with locking mechanism

## Package Includes

PCN Catheter with metal straightner, metal trocar, urine bag connector & scalpel blade

## PCN CATHETER WITH NEEDLE

### Size (in Fr)

6, 7, 8, 8.5, 9, 10,  
12, 14, 16 & 18

### Length (in cm)

20, 22, 25, 30, 35  
& 40

## PCN CATHETER WITH NEEDLE & LOCKING MECHANISM

### Size (in Fr)

6, 7, 8, 10, 12, 14 &  
16

### Length (in cm)

20, 22, 25, 30, 35  
& 40

