

PCN SET - PIGTAIL WITH LOCKING MECHANISM (NEPHROSTOMY CATHETER SET)

Used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage.

Advantages

Kit provides entire devices for performing Nephrostomy procedure

Set Includes

- IP Needle
- Guidewire J tip PTFE coated
- Fascial Dilators
- PCN Catheter with Locking Mechanism
- Urine bag connector
- Scalpel blade
- Catheter fixing disc

Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

Length (in cm)

20, 22, 25, 30, 35 & 40

