


## PTA Balloon Dilator

### Intended Use:

Used for transluminal dilation



### Product Specifications

- Shaft size: 5 Fr
- Catheter length : 80cm
- Balloon length : 4 cm
- Balloon expansion diameter : 6, 8 & 10 mm
- Nominal pressure : 7 atm
- Rate burst pressure : 14 atm
- Use guidewire : 0.035 inch

### Key Features and Benefits:

- Controlled inflation pressure allows dilating force at a stricture or narrowing
- Radiopaque marker bands enhance visualization of balloon position under fluoroscopy