


# PTA BALLOON DILATOR

Used for transluminal dilation.

## Key Features

- Controlled inflation pressure allows dilating force at a stricture or narrowing
- Radiopaque marker bands enhance visualization of balloon position under fluoroscopy



## Product Specifications

- Shaft Size: 5 Fr
- Catheter Length: 80 cm
- Balloon Length: 4 cm
- Balloon Expansion Diameter: 6, 8 & 10 mm
- Nominal Pressure: 7 atm
- Rate Burst Pressure: 14 atm
- Use Guidewire: 0.035 inch