

DENSSIFLO[®]

PTBD CATHETER SET - MICRO PUNCTURE

Used for percutaneous placement of a catheter in bile duct upto the duodenum for biliary drainage.

Set Includes

- Scalpel blade
- Chiba needle 21 G
- 0.018 inch Cope mandril wire
- 6 Fr coaxial introducer sheath
- PTFE guidewire 0.035 inch J Tip
- Fascial dilators
- Metal & plastic straightner
- PTBD catheter with metal band & hydrophilic coating Catheter fixing disc
- Drainage bag connector



Size (in Fr)

6, 7, 8, 8.5, 9, 10, 12, 14 & 16

Length (in cm)

30, 35, 40, 45 & 50

Advantages

Kit provides entire devices for performing PTBD procedure with minimal 21 G Puncture