

# DENSSIFLO<sup>®</sup>

## PTBD CATHETER WITH METAL BAND

Used for percutaneous placement of a trans hepatic biliary drainage catheter upto the duodenum for Internal - external drainage of bile.

### Key Features

- The pigtail tip configuration helps in the retention of the distal tip
- Eye holes for better drainage
- Hydrophilic coating for easy insertion
- Available with 16 or 32 Holes
- Metal band present near to last distal hole helps in correct placement of the catheter



### Type

Pigtail

### Size (in Fr)

6, 7, 8, 8.5, 9, 10, 12, 14 & 16

### Length (in cm)

30, 35, 40, 45 & 50

### Package Includes

PTBD catheter pigtail with metal band and hydrophilic coating, metal & plastic straightener and drainage bag connector.