


Pediatric Biliary “T” Tube

Intended Use:

Used for insertion into the CBD during complicated hepatobiliary surgery for external drainage.



Size (in Fr)
6, 8 & 10

Length (in cm)
15

Key Features and Benefits:

- Soft material for patient comfort
- Atraumatic design for use in infants and children
- Soft T section allows easy placement and removal of the device
- Flop C cut channels avoid coagulation & provides excellent drainage
- Optimal radiopacity for better monitoring of the position
- Syringe adaptor for post-operative dye study