


Pediatric Suprapubic Catheter Malecot

Intended Use:

Used to be inserted percutaneously into the urinary bladder and perform external drainage of urine from the urinary bladder.



Size (in Fr)
8, 10 & 12

Length (in cm)
15

Key Features and Benefits:

- Integral flower design matches with the catheter size when closed and cause less trauma during insertion
- Malecot flower to avoid migration and enhanced drainage
- Suited for small bladders as this occupies less space as compared to balloon or Pigtail catheter
- Trocar & catheter designed to prevent buckling of catheter during the procedure
- Smooth surface for ease of introduction
- Radiopaque for radiologic confirmation of the position