

PCN Set - Malecot

Intended Use:

Used for percutaneous placement of a pigtail catheter in the renal pelvis to perform external drainage of urine.



Set Includes

- IP Needle
- Guidewire J tip PTFE coated
- Fascial Dilators
- Malecot catheter
- Urine bag connector
- Scalpel blade
- Catheter fixing disc

Size (in Fr)

8, 10, 12, 14, 16, 18, 20 & 22

Length (in cm)

20, 25, 30 & 35

Benefit:

Kit provides entire devices for performing Nephrostomy procedure