


PrimacoreTM Biopsy Gun - Semi Automatic

Designed for **smooth** and **efficient** collection of **exact, high quality** core samples from soft tissues and organs.



SMOOTH access:

- Ergonomic handle with multiple grip options
- Excellent tip geometry to reduce penetration force and tissue displacement
- Exemplary tip visibility under ultrasound guidance

EFFICIENT operation:

- Convenient adjustable throw lengths – 10mm and 20mm
- Smooth stylet insertion
- Siliconized needles to reduce needle deflection
- High Stylet stiffness avoids needle bending

EXACT, High-quality sample collection:

- Penetration depth indicator – confirms primed penetration depth
- Semi-automatic firing – aids visual verification of sample notch placement
- Superlative tactile feedback of tissue penetration and collection

Set Includes:

- Primacore™ Biopsy Instrument
- Compatible Coaxial
- Trocar and Blunt Tip Stylet
- Depth Stopper
- 10 mm Adapter

Product Specifications

Size (G)	Length (cm)
14	10, 16 & 20
16	10, 16 & 20
18	10, 16, 20 & 25
20	10, 16 & 20


Clinicians Trust. Patients Wellbeing



Blue Neem Medical Devices Private Limited
Plot Nos 270 & 271, Road No 5
Harohalli Industrial Area, Phase II
Kanakapura Taluk, Ramanagara Dist - 562112
Karnataka, India



+91 97399 72855



enquiry@blueneem.com



www.blueneem.com

360° Product
Experience

