


# Primacore™

Biopsy Gun - Semi Automatic

Designed for **smooth** and **efficient** collection of **exact, high quality** core samples from soft tissues and organs.



**Smooth Stylet Insertion**  
Ensures minimal or no tissue displacement



**Echogenic Needle Tip**  
Excellent tip visibility under ultrasound guidance

**Depth Markings**  
To facilitate the precise placement of the needle at the desired depth

**Penetration Depth Indicator**  
Confirms primed penetration depth of 10mm or 20mm

**Ergonomic Handle Design**  
Provides multiple grip options for easy operation

**Fire-ready Indicator**  
Designed to reduce the risk of premature firing

## SMOOTH access:

- Ergonomic handle with multiple grip options
- Excellent tip geometry to reduce penetration force and tissue displacement
- Exemplary tip visibility under ultrasound guidance

## EFFICIENT operation:

- Convenient adjustable throw lengths – 10mm and 20mm
- Smooth stylet insertion
- Siliconized needles to reduce needle deflection
- High Stylet stiffness avoids needle bending

## EXACT, High-quality sample collection:

- Penetration depth indicator – confirms primed penetration depth
- Semi-automatic firing – aids visual verification of sample notch placement
- Superlative tactile feedback of tissue penetration and collection

## Set Includes:

- Primacore™ Biopsy Instrument
- Compatible Coaxial
- Trocar and Blunt Tip Stylet
- Depth Stopper
- 10 mm Adapter

## Product Specifications

Size (G)	Length (cm)
14	10, 16 & 20
16	10, 16 & 20
18	10, 16, 20 & 25
20	10, 16 & 20

**BLUENEEM**®  
CLINICIANS' TRUSTED PARTNER



**Blue Neem Medical Devices Private Limited**  
Plot Nos 270 & 271, Road No 5  
Harohalli Industrial Area, Phase II  
Kanakapura Taluk, Ramanagara Dist - 562112  
Karnataka, India



+91 97399 72855



enquiry@blueneem.com



www.blueneem.com

360° Product Experience

