



**PRIME WIRE<sup>®</sup>**  
PTFE NITINOL GUIDEWIRE WITH HYDROPHILIC TIP



**Type**  
Straight Tip

**Size (in inch)**  
0.025, 0.032, 0.035  
& 0.038

**Length (in cm)**  
150

### Intended Use:

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures.

### Key Features:

- Designed to combine the benefits of multiple specialty guidewires in one guidewire for urological procedures
- 5 cm hydrophilic floppy tip provides easy passage in tortuous anatomy
- Teflon coating provides low friction navigation and easy advancement of devices

### Advantages:

- One wire for all urology procedures
- Access & stenting with single wire
- Avoid exchange of wires
- Less procedural time
- Easy access in tortuous anatomy