


Prime Wire[®]

PTFE Nitinol Guidewire

Intended Use:

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures



Types
Straight & Angled Tip

Size (in inch)
0.032, 0.035 & 0.038

Length (in cm)
150 & 260

Key Features and Benefits:

- **Hybrid guidewire:** One wire for all diagnostic or interventional procedures, simplifies procedure and reduces inventory
- **Nitinol core with Teflon coating:** Provides ideal stiffness with low friction, making it easy to access and deploy devices with a single wire
- **Hydrophilic Floppy tip (5cm):** Easy access in tortuous anatomy.
- **Single Guidewire:** Avoid exchange of wires, less procedural time and reduced risk of complications