





PRIME WIRE™

PTFE NITINOL GUIDEWIRE WITH HYDROPHILIC TIP

Straight Tip
Angled Tip



ADVANTAGE OF PRIME WIRE™

Each urological procedure has its unique need. The Prime Wire™ has balanced construction to be the primary guidewire for all major urological procedures like URS, RIRS and PCNL

PRODUCT SPECIFICATIONS

| Product Code | Size (inch) | Length (cm) | Type |
|--------------|-------------|-------------|----------|
| BPTN032150HS | 0.032 | 150 | Straight |
| BPTN035150HS | 0.035 | 150 | Straight |
| BPTN038150HS | 0.038 | 150 | Straight |
| BPTN032150HA | 0.032 | 150 | Angled |
| BPTN035150HA | 0.035 | 150 | Angled |
| BPTN038150HA | 0.038 | 150 | Angled |

BLUENEEM®



Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5
Harohalli Industrial Area, Phase II
Kanakapura Taluk, Ramanagara Dist-562112
Karnataka, India
www.bluneem.com



T: +91 80 2976 1335 / 36




M: +91 97399 72855



E: marketing@bluneem.com





SINGLE WIRE FOR ALL YOUR PROCEDURAL NEEDS

PRIME WIRE[™] IMPROVES EFFICIENCY AND SAVES TIME BY REDUCING STEPS INVOLVED WITH EXCHANGE OF WIRES

ACCESS

➤ Prime Wire[™] provides a flexible and slippery tip to pass the obstruction and negotiate tortuous anatomy

MAINTAIN ACCESS

➤ Its coil construction provides sufficient retention friction making it an ideal safe wire to maintain access

DEVICE DEPLOYMENT

➤ Stiff shaft with PTFE coating for easy device manipulation and placement without kinking