

## Re-entry Malecot Drainage Catheter



### Package Includes

Re-entry Malecot Catheter with plastic straightener & urine bag connector

### Size (in Fr)

10/6, 12/6, 14/6, 16/8 & 18/8

### Length (in cm)

45

### Intended Use:

Used to be inserted percutaneously into the renal pelvis upto the ureter to perform external drainage of urine from the renal pelvis and internal drainage into the bladder.

### Key Features:

- The Malecot flower design is suitable for drainage of thick viscous fluids
- Integral flower design matches with the catheter size when closed & cause less trauma during insertion
- Used as a nephrostomy catheter and also as ureteral stent