

RE-ENTRY MALECOT DRAINAGE CATHETER

Used to be inserted percutaneously into the renal pelvis upto the ureter to perform external drainage of urine from the renal pelvis and internal drainage into the bladder.

Key Features

- The Malecot flower design is ideally suitable for the drainage of thick viscous fluids
- Integral flower design matches with the catheter size when closed and cause less trauma during insertion
- Used as a nephrostomy catheter and also as ureteral stent

Package Includes

Re-entry Malecot Catheter with plastic straightener & urine bag connector

Size (in Fr)

8/4, 10/4, 10/5, 10/6, 12/6, 14/6, 14/7, 16/6, 16/7, 16/8 & 18/8

Length (in cm)

30, 40 & 45

