

# S QUADRA FLATWIRE

## (STAINLESS STEEL STONE BASKET HELICAL WITH TIP)

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

### Key Features

- The flat wire configuration provides equal strength and securing
- Is designed to securely grasp and manipulate calculi in the ureter as well as in the renal pelvis
- Adjustable handle knob allows controlled and assured basket expansion
- In case of stone impregnation in the basket during the procedure, ergonomic handle

### Size (in Fr)

2.5, 3, 4 & 5

### Length (in cm)

70, 90, 110 & 120

