

# S QUADRA HELICAL

## (STAINLESS STEEL STONE BASKET HELICAL WITH TIP)

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

### Key Features

- Helical design helps in capturing small stones and fragments effectively
- Adjustable handle knob allows controlled and assured basket expansion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly
- Atraumatic tip

### Size (in Fr)

2.5, 3, 4 & 5

### Length (in cm)

70, 90, 110 & 120

