

KIDSOFT SIGMA STENT (URETERAL STENT SET)



Set Includes

Guidewire, Stent,
Dilator & Connector
Hub

Size (in Fr)

2.8, 3, 3.5 & 4

Length (in cm)

10/15

Intended Use:

Intended to drain urine after pyeloplasty procedure

Key Features:

- The central 'pigtail' remains in the pelvis while the narrow distal extensions stents the anastomosis and the ureter. The distal ureteric end can be trimmed if needed
- 90 degree angle avoids stent migration
- The proximal extension tubing can be used for external drainage of urine from the kidney or can be blocked by using the clamp provided
- Prevents post-operative bladder spasms
- Can be removed on an out-patient basis
- Avoids subsequent surgical intervention for stent removal with its attendant costs and potential for infections
- Extension dia 6Fr