


## Single Step<sup>®</sup> Ureteral Dilator



### Size (in Fr)

8/6, 10/5, 10/6, 10/8, 12/5  
12/6, 12/7, 12/8, 14/5, 14/6  
14/7, 14/8, 16/6 & 16/7

### Length (in cm)

70

### Intended Use:

Used for one step dilation of the ureter through natural orifice prior to a diagnostic or interventional procedure.

### Key Features:

- Long tapered tip avoids use of sequential dilation
- Smooth surface for atraumatic entry during dilation
- Also available with hydrophilic coating for ease of introduction
- Optimal radiopacity for better visualization