

# SINGLE STEP<sup>®</sup> FASCIAL DILATOR

## (LONG TAPER FASCIAL DILATOR)

Used for single step fascial dilation, intended for stretching or enlarging the fascial tissue covering a cavity, tract or opening prior to a diagnostic or interventional procedure.

### Key Features

- Long taper tip for ease of introduction & avoids use of sequential dilation
- Smooth surface for atraumatic entry during dilation
- Optimal radiopacity for better visualization
- Also available with hydrophilic coating for ease of introduction



### Size (in Fr)

6, 8, 10, 12, 14 & 16

### Length (in cm)

20 & 22