

SINGLE STEP[®] URETERAL DILATOR

Used for one step dilation of the ureter through natural orifice prior to a diagnostic or interventional procedure.

Key Features

- Long tapered tip avoids use of sequential dilation
- Smooth surface for atraumatic entry during dilation
- Also available with hydrophilic coating for ease of introduction
- Optimal radiopacity for better visualization

Size (in Fr)

8/6, 10/5, 10/6, 10/8, 12/5, 12/6, 12/7
12/8, 14/5, 14/6, 14/7, 14/8, 16/6 & 16/7

Length (in cm)

70

