

## TURP Electrodes

### Intended Use:

The above products are manufactured compatible to resectoscopes of Karl Storz, Wolf, Olympus and ACMI

### Types:

#### Cutting Loop

Used during transurethral resection of the prostate (TURP) - Single / Double stem in 24/27 Ch

#### Collins Knife

Used for ureteral meatotomy transurethral incision of the prostate (TUIP) - Single / Double stem

#### Cold Knife

For visual internal urethrotomy for strictures. Single or double stem, straight, wavered, half moon and sickle

#### Bipolar Loops

For Storz and Olympus. Single/ Double stem

#### Grooved Roller Electrode, Plain Roller Electrode and Ball Electrode


Used for coagulation - Single/ Double Stem, 3 & 5 mm

#### Monopolar Flexible Bugbee Electrodes


Used for coagulation

#### Bipolar Buttons

For Storz and Olympus. Single/ Double Stem




**Cutting Loop Standard**




**Collins Knife**




**Cold Knife - Straight**



**Ball Electrode**




**Cold Knife - Half Moon**



**Bipolar Loops**



**Plain Roller**



**Bugbee Electrodes**