

TURP ELECTRODES

Used for cutting or coagulation during the TURP Procedure.

Key Features

The above products will be manufactured compatible to resectoscopes of Karl Storz, Wolf, Olympus & ACMI

Key Features

Cutting Loop

Used during transurethral resection of the prostate (TURP) - Single / Double stem in 24/27 Ch

Collins Knife

Used for ureteral meatotomy transurethral incision of the prostate (TUIP) - Single / Double stem

Cold Knife

For visual internal urethrotomy for strictures. Single or double stem, straight, wavered, half moon & sickle

Bipolar Loops

For Storz & Olympus. Single/ Double stem

Grooved Roller Electrode, Plain Roller Electrode & Ball Electrode

Used for coagulation - Single/ Double Stem, 3 & 5 mm

Monopolar Flexible Bugbee Electrodes


Used for coagulation

Bipolar Buttons


For Storz & Olympus. Single/ Double Stem




Cutting Loop Standard




Collins Knife




Cold Knife - Straight




Ball Electrode



Cold Knife - Half Moon



Bipolar Loops



Plain Roller