

TRACER[™]

HYDROPHILIC GUIDEWIRE



Types

- Straight Tip
- Angled Tip

Size (in inch)

0.025, 0.032, 0.035
& 0.038

Length (in cm)

150

Intended Use:

Used to gain smooth and easy access to desired anatomical location and to assist in placement of medical devices.

Key Features:

- Hydrophilic coating provides easy passage in tortuous anatomy and through lesions
- Radiopaque polymer jacket enhances fluoroscopic visualization
- Kink resistant nitinol core for optimal performance