


REACH WITH SURETY

Hydrophilic guidewire, designed to reach the desired anatomical location and support device deployment

ADVANTAGES OF TRACER™

- Hydrophilic coating provides excellent lubricity and smooth navigation
- 1 : 1 torque ratio from core to tip providing optimal feedback and controlled navigation
- Kink resistant Nitinol core provides optimal strength for easier navigation through tortuous anatomy and device deployment
- Radiopacity of the wire provide excellent visibility



| PRODUCT SPECIFICATIONS | | | |
|------------------------|-----------|-------------|----------|
| Product Code | Size (in) | Length (cm) | Type |
| BHWN018150ST | 0.018 | 150 | Straight |
| BHWN025150ST | 0.025 | 150 | Straight |
| BHWN032150ST | 0.032 | 150 | Straight |
| BHWN035150ST | 0.035 | 150 | Straight |
| BHWN038150ST | 0.038 | 150 | Straight |
| BHWN018150AT | 0.018 | 150 | Angled |
| BHWN025150AT | 0.025 | 150 | Angled |
| BHWN032150AT | 0.032 | 150 | Angled |
| BHWN035150AT | 0.035 | 150 | Angled |
| BHWN038150AT | 0.038 | 150 | Angled |

BNMD/MKG/LFT/BN/73 | Rev.01 | 20-May-2021

BLUENEEM







M: +91 97399 72855



E: marketing@blueneem.com



Blue Neem Medical Devices Pvt Ltd Plot Nos 270 & 271, Road No 5 Harohalli Industrial Area, Phase II Kanakapura Taluk, Ramanagara Dist - 562112 Karnataka, India