

Tracer®

Hydrophilic Guidewire

Intended Use:

Used to gain smooth and easy access to desired anatomical location and to assist in placement medical devices

**Types**

Straight & Angled Tip

Size (in inch)

0.025, 0.032, 0.035 & 0.038

Length (in cm)

150, 180 & 260

Key Features and Benefits:

- Hydrophilic coating provides easy passage in tortuous anatomy and through lesions
- Radiopaque polymer jacket enhances fluoroscopic visualization
- Kink resistant nitinol core for optimal performance