




TRACERTM

HYDROPHILIC GUIDEWIRE

REACH WITH SURETY

Hydrophilic guidewire, designed to reach the desired anatomical location and support device deployment

ADVANTAGES OF TRACERTM

- Hydrophilic coating provides excellent lubricity and smooth navigation
 - 1 : 1 torque ratio from core to tip providing optimal feedback and controlled navigation
 - Kink resistant Nitinol core provides optimal strength for easier navigation through tortuous anatomy and device deployment
 - Radiopacity of the wire provide excellent visibility
- 

KEY FEATURES



1 Radiopaque wire and tip for fluoroscopic visibility and controlled navigation under fluoroscopic guidance

2 Flexible & tapered tip for atraumatic and controlled navigation

3 Hydrophilic coating for low friction, steerability and trackability

4 Nitinol core - kink resistant for easy device deployment

PRODUCT SPECIFICATIONS

Product Code	Size (in)	Length (cm)	Type
BHWN025150ST	0.025	150	Straight
BHWN032150ST	0.032	150	Straight
BHWN035150ST	0.035	150	Straight
BHWN038150ST	0.038	150	Straight
BHWN025150AT	0.025	150	Angled
BHWN032150AT	0.032	150	Angled
BHWN035150AT	0.035	150	Angled
BHWN038150AT	0.038	150	Angled

BNMD/MKG/LFT/BN/73 | Rev.01 | 20-May-2021

BLUENEEM[®]

CE 0123



Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5
Harohalli Industrial Area, Phase II
Kanakapura Taluk, Ramanagara Dist - 562112
Karnataka, India



T: +91 80 2976 1335 / 36



M: +91 97399 72855



E: marketing@blueneem.com