


## Triflex Basket

### Intended Use:

Used in an endoscopic procedure for removal of urinary stones or foreign objects from calyx



**Size (in Fr)**  
2.2 & 2.4

**Length (in cm)**  
115 & 120

### Key Features and Benefits:

- Flexible design shaft which can angulate up to 270° as per the flexible scope direction during RIRS procedure
- Kink resistant braided shaft provides durability, making removal of several fragment possible
- Triflex design holds like a basket and releases like a grasper
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- Useful for stone repositioning during RIRS procedures