

Triflex Basket



Size (in Fr)
2.2 & 2.4

Length (in cm)
115 & 120

Intended Use:

Triflex Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from calyx.

Key Features:

- Flexible design shaft which can angulate up to 270 degrees as per the flexible scope direction during RIRS procedure
- Kink resistant braided shaft provides durability, making removal of several fragments possible
- Triflex design holds like a basket and releases like a grasper
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- Useful for stone repositioning during RIRS Procedures