


Ureteral Balloon Dilator



Product Specifications

- Shaft Size: 5 Fr
- Balloon Diameter: 5 mm
- Balloon Length: 4 cm
- Catheter Length: 75 cm
- RBP: 14 atm
- Use Guidewire: 0.035 inch

Intended Use:

Used for dilation of narrow ureteral segments.

Key Features:

- Controlled inflation pressure allowing for dilating force at a stricture or ureteral narrowing
- Radiopaque marker bands enhance visualization of balloon position under fluoroscopy