


KIDSOFT URETERAL STENT



Package Includes
Stent with pusher

Types
Both Ends Open

Size (in Fr)
2.8, 3, 3.5 & 4

Length (in cm)
10, 12, 14 & 16

Intended Use:

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder.

Key Features:

- Polyurethane material enables smooth surface & low encrustation tendency
- Pigtail Loop avoids slippage & has very good retention
- Radiopaque for better visualization
- Indwelling period is 30 days

Also Available In Set Which Includes

Stent with Pusher & Guidewire