


Relief CIC[®]

Hydrophilic Intermittent Catheter

Intended Use:

Used to drain urine from bladder



Key Features and Benefits:

- Hydrophilic coated for easy insertion and removal, without using Lubricating Gel
- Highly lubricious, sustained hydrophilic coating
- Clear catheter tube with balanced flexibility and stiffness
- Eye holes for better drainage

Type	Size (in Fr)	Length (in cm)
Male	10, 12, 14, 16 & 18	35
Female	10, 12, 14, 16 & 18	15