


**Suction.
Evacuation.
Collection.**

WortX™ facilitates ureteral dilation, ensuring smooth and efficient scope access to the ureter. It enhances safety and effectiveness during RIRS procedures, contributing to an improved stone-free rate (SFR).



**Flexible sheath
for easy access**



**Length markings
for accurate
placement**



**Hydrophilic coating
for easy insertion**



**Tapered tip for
smoother insertion**



Key Features And Benefits



Flexible Sheath

10cm flexible tip for easy kidney access with a scope,
a rigid shaft for scope support



Suction Vent

Integrated suction continuously removes fluids and debris,
enhancing visibility and shortening procedure time




Air-Tight Seal

The bottle features an airtight seal to ensure
effective suction throughout the procedure



Graduated Markings

Markings on the side measure the
volume of the fluid



Collector

Collects debris and prevents it from entering
the suction machine

Set Includes

1. OD of the UAS – 12 Fr
2. ID of the UAS – 10 Fr
3. Length of the UAS – 35, 40, 45 & 50 cm
4. Use Guidewire Size – 0.035 Inch
5. Stone Collection Bottle – 200 ml
6. Silicon Tube – Secure Fit for
7. Stone Collection Bottle and Device

BLUENEEM®
Clinicians Trust. Patients Wellbeing



Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5
Harohalli Industrial Area, Phase II
Kanakapura Taluk
Ramanagara Dist - 562112
Karnataka, India



+91 97399 72855



enquiry@blueneem.com