

Accura Kit - Percutaneous Access Set

Intended Use:

Used for percutaneous access with 21G needle and subsequently introducing a 0.038 inch guidewire, offering a minimally invasive entry point for Interventional Radiology procedures



Set Includes

- 6Fr introducer sheath
- 21G x 15cm chiba needle with echogenic tip
- 0.018 inch nitinol cope mandril wire

Key Features and Benefits:

- A radiopaque band located at the tip of the sheath helps to identify precise location of the sheath's distal tip for accurate positioning
- The hydrophilic coating on sheath when activated becomes highly lubricious, enabling smooth insertion by reducing friction
- Sheath lumen allows both 0.018 inch cope mandril wire and 0.038 inch wire together
- Set includes a chiba needle with echogenic tip which enhances visibility under ultrasound