


Amplatz Renal Dilation Set

Intended Use:

Intended for stretching or enlarging the fascial tissue prior to a nephrostomy procedure and to introduce the amplatz sheath which is intended for percutaneous dilation of tissue surrounding a dilator previously placed and used as a working channel during urological procedures.



Set Includes:

- 8 Fr Guiding Catheter
- 10 Fr Radiopaque Dilation Sheath
- Radiopaque Fascial dilators from 6 Fr – 10 Fr of length 20 cm
- Radiopaque Renal dilators from 12 Fr – 30 Fr of length 30 cm
- Teflon Amplatz sheath for dilators 24, 26, 28 & 30 Fr of length 17 cm
- Packed in tray