

# Amplatz Renal Dilator (Renal Dilator With Amplatz Sheath)

**Intended Use:**

Used for dilation of a percutaneous tract and setup of a working channel into the kidney for PCNL procedures by placing a suitable amplatz sheath



**Key Features and Benefits:**

- Gradual taper designed for atraumatic placement and maneuverability
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization

Type	Size(OD of Dilator in Fr)	Length (in cm)
Dilator	12, 14, 16, 18, 20, 22, 24 26, 28 & 30	30

Type	Size(ID of Sheath in Fr)	Length (in cm)
Amplatz Sheath	12, 14, 16, 18, 20, 22, 24 26, 28 & 30	17