

## Amplatz Sheath

### Intended Use:

Intended for percutaneous dilation of tissue surrounding a dilator previously placed. It is used as a working channel to maintain a previously established nephrostomy tract for introduction of instruments and catheters during urological procedures



#### Size (in Fr)

12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

#### Length (in cm)

17 & 20

#### Key Features and Benefits:

- Optimal bevel tip for atraumatic entry
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization