

Amplatz Ultra Stiff® Guidewire

Intended Use:

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures



Type	Size (in inch)	Length (in cm)
J Tip	0.035	80

Key Features and Benefits:

- PTFE coated stainless steel mandril with 5 cm atraumatic floppy tip for less trauma
- Smooth transition from stiff mandril to the flexible tip
- The stiffness of the mandril design provides support for the advancement of devices such as dilators, catheters and instruments