


## Dalela Triple J Ureteral Stent

### Intended Use:

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder



### Size (in Fr)

3, 3.5, 3.7, 3.8, 4, 4.5, 4.7, 4.8, 5, 5.5, 6, 6.5, 7 & 8

### Length (in cm)

22, 24 & 26

### Key Features and Benefits:

- **Advanced 'J' Bend Design:** Minimises irritation by reducing contact with the sensitive trigone
- **Thermosensitive Advanced PU Material:** Ensures easy deployment and enhanced comfort at body temperature
- **Reduction of Dysuria and Stent-Induced Symptoms:** Provides greater comfort for patients
- **Long-Term Indwelling Design with Radiopaque Material:** Enables better visualization during imaging
- **Easy Removal with Visible Lower Coil:** Simplifies the stent removal process