


Intended Use:

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder



Types

Both Ends Open
One End Closed
& With or Without String

Size (in Fr)

5, 5.5, 6, 6.5, 7 & 8

Length (in cm)

12, 14, 16, 18, 20, 22, 24
26, 28 & 30

Package Includes
Stent with pusher

Also available in set, which Includes
Stent with pusher and guidewire

Key Features and Benefits:

- Silicone blend polyurethane material for enhanced patient comfort
- Optimal loop shape and high memory retention reduces the risk of migration
- Optimal radiopacity for better visualization