

Initial Puncture Needle - Two Part

Intended Use:

Intended to make the initial percutaneous puncture to access the kidney and to enable the introduction of a guidewire during urological procedures (PCN, PCNL, Mini PERC and Micro PERC)



Type	Size (in G)	Length (in cm)
Two part	14, 16 & 18	10, 15, 20, 23 & 25

Key Features and Benefits:

- **Siliconized needle:** Reduces trauma for a smoother procedure
- **Echogenic tip:** Enhances visibility under ultrasound
- **1cm shaft marking:** Provides precise depth control
- **Locking hub:** Ensures stability for easier initial puncture during the procedure. Locking hub helps in holding the inner needle with other cannula together, which helps in easy initial puncture